


Assembly Instructions


ELEVATE

Electric Lift Chair




WARNING


FAILURE TO FOLLOW THESE BASIC STEPS CAN LEAD TO DAMAGE THAT IS NOT COVERED BY YOUR WARRANTY



Read Instructions carefully. Check no parts are missing.




Lay out all of the components on a soft surface as if you were preparing to do a jigsaw puzzle.




Carefully identify each component, especially those that are similar. The most common mistake is getting things back to front.



Separate and identify each component from the hardware pack and rest each one in its proper location before fastening.




DO NOT fully tighten screws (leave slightly loose). When the item is fully assembled tighten the screws - but do not over tighten.



If glue is required **DO NOT** set until you are certain that the parts belong together.

- A. BACK CUSHION 1
- B. SEAT CUSHION 1
- C. VELCRO 1

CAUTION: NEED AT LEAST TWO PEOPLE TO COMPLETE THE ASSEMBLY.



1. CONNECT THE BACK CUSHION(A) TO THE SEAT (B) BY INSERTING CLIP.
2. PULL DOWN THE VELCRO (C) ON THE BACK COVER AND ATTACH IT UNDER THE BACK OF THE CHAIR.